

Bovine Troponin I Type 3, Cardiac ELISA Kit | cTnI elisa kit

Bovine Troponin I Type 3, Cardiac ELISA Kit

CATALOG #
AAA90038

REACTIVITY
Bovine


PRODUCT SYNONYMS

Troponin I Type 3; Cardiac; Bovine Troponin I Type 3; Cardiac ELISA Kit; cTnI elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	1.0 ng/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for cTnI elisa kit

For samples: Serum, plasma, cell culture supernatants, body fluid and tissue homogenate

INTENDED USE This CTNI ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Bovine CTNI. This ELISA kit for research use only!

PRINCIPLE OF THE ASSAY CTNI ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-CTNI antibody and an CTNI-HRP conjugate. The assay sample and buffer are incubated together with CTNI-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the CTNI concentration since CTNI from samples and CTNI-HRP conjugate compete for the anti-CTNI antibody binding site. Since the number of sites is limited, as more sites are occupied by CTNI from the sample, fewer sites are left to bind CTNI-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The CTNI concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/ctni-bovine-ovine-elisa-kits/troponin-i-type-3-cardiac/90038>