

Goat Insulin Autoantibody ELISA Kit | IAA elisa kit

Goat Insulin Autoantibody ELISA Kit

CATALOG #
AAA90563

REACTIVITY
Goat


PRODUCT SYNONYMS

Insulin Autoantibody; Goat Insulin Autoantibody ELISA Kit; IAA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Goat
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Sandwich
Sensitivity	1.0 ng/mL.
Intended Uses	This IAA ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Goat IAA. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IAA elisa kit

Principle of the Assay: IAA ELISA kit applies the quantitative sandwich enzyme immunoassay technique. The microtiter plate has been pre-coated with a monoclonal antibody specific for IAA. Standards or samples are then added to the microtiter plate wells and IAA if present, will bind to the antibody pre-coated wells. In order to quantitatively determine the amount of IAA present in the sample, a standardized preparation of horseradish peroxidase (HRP) -conjugated polyclonal antibody, specific for IAA are added to each well to "sandwich" the IAA immobilized on the plate. The microtiter plate undergoes incubation, and then the wells are thoroughly washed to remove all unbound components. Next, substrate solutions are added to each well. The enzyme (HRP) and substrate are allowed to react over a short incubation period. Only those wells that contain IAA and enzyme-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The IAA concentration in each sample is interpolated from this standard curve.

Product Categories/Family for IAA elisa kit

Immunology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/iaa-goat-elisa-kits/insulin-autoantibody/90563>