

Canine Salusin Beta ELISA Kit | Salusin B elisa kit

Canine Salusin Beta ELISA Kit

CATALOG #
AAA90570

REACTIVITY
Canine


PRODUCT SYNONYMS

Salusin Beta; Canine Salusin Beta ELISA Kit; Salusin B elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Sandwich
Detection Range	100-2500pg/mL
Sensitivity	1.0pg/mL
Intended Uses	This Salusin Beta ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine Salusin Beta. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Salusin B elisa kit

Principle of the Assay: Salusin Beta ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-Salusin Beta antibody and an Salusin Beta-HRP conjugate. The assay sample and buffer are incubated together with Salusin Beta-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the Salusin Beta concentration since Salusin Beta from samples and Salusin Beta-HRP conjugate compete for the anti-Salusin Beta antibody binding site. Since the number of sites is limited, as more sites are occupied by Salusin Beta from the sample, fewer sites are left to bind Salusin Beta-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The Salusin Beta concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/salusin-b-canine-elisa-kits/salusin-beta/90570>