

Bovine Pro Atrial Natriuretic Peptide ELISA Kit | Pro-ANP elisa kit

Bovine Pro Atrial Natriuretic Peptide ELISA Kit

CATALOG # AAA90765	REACTIVITY Bovine
------------------------------	-----------------------------


PRODUCT SYNONYMS

Pro Atrial Natriuretic Peptide; Bovine Pro Atrial Natriuretic Peptide ELISA Kit; Pro-ANP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	1.0 ng/mL.
Intended Uses	This Pro-ANP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Bovine Pro-ANP. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Pro-ANP elisa kit	Principle of the Assay: Pro-ANP ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-Pro-ANP antibody and an Pro-ANP -HRP conjugate. The assay sample and buffer are incubated together with Pro-ANP -HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the Pro-ANP concentration since Pro-ANP from samples and Pro-ANP -HRP conjugate compete for the anti-Pro-ANP antibody binding site. Since the number of sites is limited, as more sites are occupied by Pro-ANP from the sample, fewer sites are left to bind Pro-ANP -HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The Pro-ANP concentration in each sample is interpolated from this standard curve.
Product Categories/Family for Pro-ANP elisa kit	Cell Biology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/pro-anp-bovine-ovine-elisa-kits/pro-atrial-natriuretic-peptide/90765>