

Canine Complement 1 Inhibitor ELISA Kit | C1INH) elisa kit

Canine Complement 1 Inhibitor ELISA Kit

CATALOG #
AAA91376

REACTIVITY
Canine


PRODUCT SYNONYMS

Complement 1 Inhibitor; Canine Complement 1 Inhibitor ELISA Kit; C1INH) elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Sandwich
Detection Range	25-500ng/mL
Sensitivity	1.0ng/mL
Intended Uses	This C1INH ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine C1INH. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for C1INH) elisa kit

Principle of the assay: C1INH ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-C1INH antibody and an C1INH-HRP conjugate. The assay sample and buffer are incubated together with C1INH-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the C1INH concentration since C1INH from samples and C1INH-HRP conjugate compete for the anti-C1INH antibody binding site. Since the number of sites is limited, as more sites are occupied by C1INH from the sample, fewer sites are left to bind C1INH-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The C1INH concentration in each sample is interpolated from this standard curve.

Product Categories/Family for C1INH) elisa kit

Immunology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/c1inh-canine-elisa-kits/complement-1-inhibitor/91376>