

Mouse alpha-Galactosyl ELISA Kit | alpha-Gal elisa kit

Mouse alpha-Galactosyl ELISA Kit

CATALOG #
AAA91423

REACTIVITY
Mouse


PRODUCT SYNONYMS

alpha-Galactosyl; Mouse alpha-Galactosyl ELISA Kit; Mouse a-Galactosyl ELISA Kit; alpha-Gal elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	1-25ng/mL
Sensitivity	0.1ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for alpha-Gal elisa kit

Intended Uses: This alphaGAL ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Mouse alphaGAL. This ELISA kit for research use only!

Principle of the Assay||alphaGAL ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-alphaGAL antibody and an alphaGAL-HRP conjugate. The assay sample and buffer are incubated together with alphaGAL-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the alphaGAL concentration since alphaGAL from samples and alphaGAL-HRP conjugate compete for the anti-alphaGAL antibody binding site. Since the number of sites is limited, as more sites are occupied by alphaGAL from the sample, fewer sites are left to bind alphaGAL-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The alphaGAL concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

