

Sheep Vitamin D Receptor ELISA Kit | VD R elisa kit

Sheep Vitamin D Receptor ELISA Kit

CATALOG #
AAA92861

REACTIVITY
Sheep


PRODUCT SYNONYMS

Vitamin D Receptor; Sheep Vitamin D Receptor ELISA Kit; VD R elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Sheep
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	0.1 ng/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for VD R elisa kit

For samples: Serum, plasma, cell culture supernatants, body fluid and tissue homogenate

INTENDED USE This VDR ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Sheep VDR. This ELISA kit for research use only!

PRINCIPLE OF THE ASSAY VDR ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-VDR antibody and an VDR-HRP conjugate. The assay sample and buffer are incubated together with VDR-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the VDR concentration since VDR from samples and VDR-HRP conjugate compete for the anti-VDR antibody binding site. Since the number of sites is limited, as more sites are occupied by VDR from the sample, fewer sites are left to bind VDR-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The VDR concentration in each sample is interpolated from this standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4262865
NCBI Official Full Name	vitamin D receptor, partial

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/vd-r-sheep-elisa-kits/vitamin-d-receptor/92861>