

## Canine Pirin ELISA Kit | PIR elisa kit

Canine Pirin ELISA Kit

CATALOG #  
**AAA93260**

REACTIVITY  
**Canine**


### PRODUCT SYNONYMS

Pirin; Canine Pirin ELISA Kit; PIR elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                                                                                                                                       |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity              | Canine                                                                                                                                                                                |
| Samples                 | Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate                                                                                                            |
| Assay Type              | Sandwich                                                                                                                                                                              |
| Detection Range         | 1.0-25ng/mL                                                                                                                                                                           |
| Sensitivity             | 0.1ng/mL                                                                                                                                                                              |
| Intended Uses           | This PIR ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine PIR. This ELISA kit for research use only, not for therapeutic applications! |
| Preparation and Storage | Store all reagents at 2-8 degree C.                                                                                                                                                   |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for PIR elisa kit

**Principle of the Assay:** PIR ELISA kit applies the quantitative sandwich enzyme immunoassay technique. The microtiter plate has been pre-coated with a monoclonal antibody specific for PIR. Standards or samples are then added to the microtiter plate wells and PIR if present, will bind to the antibody pre-coated wells. In order to quantitatively determine the amount of PIR present in the sample, a standardized preparation of horseradish peroxidase (HRP) -conjugated polyclonal antibody, specific for PIR are added to each well to "sandwich" the PIR immobilized on the plate. The microtiter plate undergoes incubation, and then the wells are thoroughly washed to remove all unbound components. Next, substrate solutions are added to each well. The enzyme (HRP) and substrate are allowed to react over a short incubation period. Only those wells that contain PIR and enzyme-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The PIR concentration in each sample is interpolated from this standard curve.

## NCBI AND UNIPROT INFORMATION

|                            |                          |
|----------------------------|--------------------------|
| NCBI GI #                  | <a href="#">49456601</a> |
| Molecular Weight           | 32,113 Da                |
| NCBI Official Full Name    | PIR, partial             |
| UniProt Protein Name       | Pirin                    |
| UniProt Gene Name          | <b>PIR</b>               |
| UniProt Synonym Gene Names | Probable quercetinase    |
| UniProt Entry Name         | PIR_HUMAN                |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/pir-canine-elisa-kits/pirin/93260>