

Porcine Aprotinin ELISA Kit | AP elisa kit

Porcine Aprotinin ELISA Kit

CATALOG #
AAA93691

REACTIVITY
Porcine


PRODUCT SYNONYMS

Aprotinin; Porcine Aprotinin ELISA Kit; AP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Detection Range	5.0-100ng/mL
Sensitivity	1.0ng/mL
Intended Uses	This AP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Porcine AP. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for AP elisa kit

Principle of the Assay: AP ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-AP antibody and an AP-HRP conjugate. The assay sample and buffer are incubated together with AP-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the AP concentration since AP from samples and AP-HRP conjugate compete for the anti-AP antibody binding site. Since the number of sites is limited, as more sites are occupied by AP from the sample, fewer sites are left to bind AP-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The AP concentration in each sample is interpolated from this standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	58005
NCBI Official Full Name	aprotinin, partial

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/ap-porcine-elisa-kits/aprotinin/93691>