

Porcine Cyclic Guanosine Monophosphate ELISA Kit | cGMP elisa kit

Porcine Cyclic Guanosine Monophosphate ELISA Kit

CATALOG #
AAA93929

REACTIVITY
Porcine


PRODUCT SYNONYMS

Cyclic Guanosine Monophosphate; Porcine Cyclic Guanosine Monophosphate ELISA Kit; cGMP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	1.0 pmol/ml.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for cGMP elisa kit

Intended Uses: This cGMP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Porcine cGMP. This ELISA kit for research use only!

Principle of the Assay: cGMP ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-cGMP antibody and an cGMP-HRP conjugate. The assay sample and buffer are incubated together with cGMP-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the cGMP concentration since cGMP from samples and cGMP-HRP conjugate compete for the anti-cGMP antibody binding site. Since the number of sites is limited, as more sites are occupied by cGMP from the sample, fewer sites are left to bind cGMP-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The cGMP concentration in each sample is interpolated from this standard curve.

Product Categories/Family for cGMP elisa kit

Signal Transduction

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/cgmp-porcine-elisa-kits/cyclic-guanosine-monophosphate/93929>