

## Canine Phospho Epidermal growth factor receptor ELISA Kit | P-EGFR elisa kit

Canine Phospho Epidermal growth factor receptor ELISA Kit

CATALOG #  
**AAA94476**

REACTIVITY  
**Canine**


### PRODUCT SYNONYMS

Phospho Epidermal growth factor receptor; Canine Phospho Epidermal growth factor receptor ELISA Kit; P-EGFR elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                                                                                                                                             |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity              | Canine                                                                                                                                                                                      |
| Samples                 | Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate                                                                                                                  |
| Sensitivity             | 0.1 ng/mL.                                                                                                                                                                                  |
| Intended Uses           | This P-EGFR ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine P-EGFR. This ELISA kit for research use only, not for therapeutic applications! |
| Preparation and Storage | Store all reagents at 2-8 degree C.                                                                                                                                                         |

### IMAGES



### ADDITIONAL INFORMATION

Product Categories/Family for P-EGFR  
elisa kit

Signal Transduction

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://stage.aaabiotech.com/p-egfr-canine-elisa-kits/phospho-epidermal-growth-factor-receptor/94476>