

Bovine Immunoglobulin G2 ELISA Kit | IgG2 elisa kit

Bovine Immunoglobulin G2 ELISA Kit

CATALOG #
AAA94526

REACTIVITY
Bovine


PRODUCT SYNONYMS

Immunoglobulin G2; Bovine Immunoglobulin G2 ELISA Kit; IgG2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Samples	Serum, Plasma, Cell Culture Supernatants, Body Fluid And Tissue Homogenate
Assay Type	Quantitative Competitive
Sensitivity	1.0 mg/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IgG2 elisa kit

Intended Uses: This IgG ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Porcine IgG. This ELISA kit for research use only!

Principle of the Assay: IgG ELISA kit applies the competitive enzyme immunoassay technique utilizing a polyclonal anti-IgG antibody and an IgG-HRP conjugate. The assay sample and buffer are incubated together with IgG-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the IgG concentration since IgG from samples and IgG-HRP conjugate compete for the anti-IgG antibody binding site. Since the number of sites is limited, as more sites are occupied by IgG from the sample, fewer sites are left to bind IgG-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The IgG concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/igg2-bovine-ovine-elisa-kits/immunoglobulin-g2/94526>