

Rat Glucose, Free ELISA Kit | GF elisa kit

Rat Glucose, Free ELISA Kit

CATALOG #
AAA94823

REACTIVITY
Rat


PRODUCT SYNONYMS

Glucose; Free; Rat Glucose; Free ELISA Kit; GF elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	0.1 ug/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GF elisa kit

Intended Uses: This FGLU ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Rat FGLU. This ELISA kit for research use only!

Principle of the Assay||FGLU ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-FGLU antibody and an FGLU-HRP conjugate. The assay sample and buffer are incubated together with FGLU-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the FGLU concentration since FGLU from samples and FGLU-HRP conjugate compete for the anti-FGLU antibody binding site. Since the number of sites is limited, as more sites are occupied by FGLU from the sample, fewer sites are left to bind FGLU-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The FGLU concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.