

Canine Collagenase, Type II ELISA Kit | Collagenase II elisa kit

Canine Collagenase, Type II ELISA Kit

CATALOG #
AAA94933

REACTIVITY
Canine


PRODUCT SYNONYMS

Collagenase; Type II; Canine Collagenase; Type II ELISA Kit; Collagenase II elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Intended Uses	This Collagenase, Type II ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine Collagenase, Type II. This ELISA kit for research use only, not for therapeutic applications
Intended Uses	This Collagenase, Type II ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine Collagenase, Type II. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Collagenase II elisa kit

Principle of the Assay: Collagenase, Type II ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-Collagenase, Type II antibody and an Collagenase, Type II-HRP conjugate. The assay sample and buffer are incubated together with Collagenase, Type II-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the Collagenase, Type II concentration since Collagenase, Type II from samples and Collagenase, Type II-HRP conjugate compete for the anti-Collagenase, Type II antibody binding site. Since the number of sites is limited, as more sites are occupied by Collagenase, Type II from the sample, fewer sites are left to bind Collagenase, Type II-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The Collagenase, Type II concentration in each sample is interpolated from this standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	203499
NCBI GeneID	171052
Molecular Weight	53,375 Da
NCBI Official Full Name	collagenase, partial
NCBI Official Synonym Full Names	matrix metalloproteinase 13
NCBI Official Symbol	Mmp13
NCBI Protein Information	collagenase 3; MMP-13; UMRCASE; collagenase-3; matrix metalloproteinase 13; matrix metalloproteinase-13
UniProt Protein Name	Collagenase 3
UniProt Gene Name	Mmp13
UniProt Synonym Gene Names	MMP-13
UniProt Entry Name	MMP13_RAT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/collagenase-ii-canine-elisa-kits/collagenase-type-ii/94933>