

Bovine Anti-Bromodeoxyuridine antibody ELISA Kit | BrdU elisa kit

Bovine Anti-Bromodeoxyuridine antibody ELISA Kit

CATALOG # AAA95353	REACTIVITY Bovine
------------------------------	-----------------------------


PRODUCT SYNONYMS

Anti-Bromodeoxyuridine antibody; Bovine Anti-Bromodeoxyuridine antibody ELISA Kit; BrdU elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	1.0 ng/mL
Intended Uses	This Anti-BrdU ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Bovine Anti-BrdU. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BrdU elisa kit

Principle of the assay: Anti-BrdU ELISA kit applies the competitive enzyme immunoassay technique utilizing BrdU antigen and an Anti-BrdU-HRP conjugate. The assay sample and buffer are incubated together with Anti-BrdU-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the Anti-BrdU concentration since Anti-BrdU from samples and Anti-BrdU-HRP conjugate compete for the BrdU antigen binding site. Since the number of sites is limited, as more sites are occupied by Anti-BrdU from the sample, fewer sites are left to bind Anti-BrdU-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The Anti-BrdU concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/brdu-bovine-ovine-elisa-kits/anti-bromodeoxyuridine-antibody/95353>