

Mouse Soluble transferrin receptor ELISA Kit | sTfR elisa kit

Mouse Soluble transferrin receptor ELISA Kit

CATALOG #
AAA95878

REACTIVITY
Mouse


PRODUCT SYNONYMS

Soluble transferrin receptor; Mouse Soluble transferrin receptor ELISA Kit; sTfR elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Quantitative Competitive
Sensitivity	1.0 ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for sTfR elisa kit

Intended Uses: This sTfER ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Mouse sTfER.

Principle of the Assay: sTfER ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-sTfER antibody and an sTfER-HRP conjugate. The assay sample and buffer are incubated together with sTfER-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the sTfER concentration since sTfER from samples and sTfER-HRP conjugate compete for the anti-sTfER antibody binding site. Since the number of sites is limited, as more sites are occupied by sTfER from the sample, fewer sites are left to bind sTfER-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The sTfER concentration in each sample is interpolated from this standard curve.

Product Categories/Family for sTfR elisa kit

[Signal Transduction](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

