

Canine Embryonic Stem MESPU30 ELISA Kit | M30 elisa kit

Canine Embryonic Stem MESPU30 ELISA Kit

CATALOG #
AAA96697

REACTIVITY
Canine


PRODUCT SYNONYMS

Embryonic Stem MESPU30; Canine Embryonic Stem MESPU30 ELISA Kit; M30 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	1.0 pg/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for M30 elisa kit

Principle of the assay: M30 ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-M30 antibody and an M30-HRP conjugate. The assay sample and buffer are incubated together with M30-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the M30 concentration since M30 from samples and M30-HRP conjugate compete for the anti-M30 antibody binding site. Since the number of sites is limited, as more sites are occupied by M30 from the sample, fewer sites are left to bind M30-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The M30 concentration in each sample is interpolated from this standard curve.

Product Categories/Family for M30 elisa kit

[Signal Transduction](#)

NCBI AND UNIPROT INFORMATION

NCBI GI #	19717569
Molecular Weight	14,749 Da
NCBI Official Full Name	M30, partial
UniProt Protein Name	M30
UniProt Gene Name	M30
UniProt Entry Name	Q8T667_ECHGR

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/m30-canine-elisa-kits/embryonic-stem-mespu30/96697>