

Porcine Intact Parathormone ELISA Kit | I-PTH elisa kit

Porcine Intact Parathormone ELISA Kit

CATALOG #
AAA96759

REACTIVITY
Porcine


PRODUCT SYNONYMS

Intact Parathormone; Porcine Intact Parathormone ELISA Kit; I-PTH elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	0.1 ng/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for I-PTH elisa kit

Principle of the assay: I-PTH ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-I-PTH antibody and an I-PTH-HRP conjugate. The assay sample and buffer are incubated together with I-PTH-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the I-PTH concentration since I-PTH from samples and I-PTH-HRP conjugate compete for the anti-I-PTH antibody binding site. Since the number of sites is limited, as more sites are occupied by I-PTH from the sample, fewer sites are left to bind I-PTH-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The I-PTH concentration in each sample is interpolated from this standard curve.

Product Categories/Family for I-PTH elisa kit

Signal Transduction

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/i-ptb-porcine-elisa-kits/intact-parathormone/96759>