

Sheep Total bile acid ELISA Kit | TBA elisa kit

Sheep Total bile acid ELISA Kit

CATALOG #
AAA96814

REACTIVITY
Sheep


PRODUCT SYNONYMS

Total bile acid; Sheep Total bile acid ELISA Kit; TBA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Sheep
Samples	Serum, Plasma, Cell Culture Supernatants, Body Fluid And Tissue Homogenate
Assay Type	Quantitative Competitive
Sensitivity	0.1 ng/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TBA elisa kit

Intended Uses: This TAU ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Rat TAU. This ELISA kit for research use only!

Principle of the Assay: TAU ELISA kit applies the competitive enzyme immunoassay technique utilizing a polyclonal anti-TAU antibody and an TAU-HRP conjugate. The assay sample and buffer are incubated together with TAU-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the TAU concentration since TAU from samples and TAU-HRP conjugate compete for the anti-TAU antibody binding site. Since the number of sites is limited, as more sites are occupied by TAU from the sample, fewer sites are left to bind TAU-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The TAU concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

