

Mouse 5-Hydroxyindoleacetic Acid ELISA Kit | 5-HIAA elisa kit

Mouse 5-Hydroxyindoleacetic Acid ELISA Kit

CATALOG #
AAA96927

REACTIVITY
Mouse


PRODUCT SYNONYMS

5-Hydroxyindoleacetic Acid; Mouse 5-Hydroxyindoleacetic Acid ELISA Kit; 5-HIAA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	0.1 mg/L.
Intended Uses	This 5HIAA ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Mouse 5HIAA. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for 5-HIAA elisa kit

Principle of the Assay: 5HIAA ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-5HIAA antibody and an 5HIAA-HRP conjugate. The assay sample and buffer are incubated together with 5HIAA-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the 5HIAA concentration since 5HIAA from samples and 5HIAA-HRP conjugate compete for the anti-5HIAA antibody binding site. Since the number of sites is limited, as more sites are occupied by 5HIAA from the sample, fewer sites are left to bind 5HIAA-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The 5HIAA concentration in each sample is interpolated from this standard curve.

Product Categories/Family for 5-HIAA elisa kit

[Epigenetics and Nuclear Signaling](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

