

Porcine Tryptophane ELISA Kit | Trp elisa kit

Porcine Tryptophane ELISA Kit

CATALOG #
AAA97241

REACTIVITY
Porcine


PRODUCT SYNONYMS

Tryptophane; Porcine Tryptophane ELISA Kit; Trp elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Detection Range	5.0-100ng/mL
Sensitivity	1.0ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Trp elisa kit

Intended Uses: This TRY ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Porcine TRY. This ELISA kit for research use only!

Principle of the Assay: TRY ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-TRY antibody and an TRY-HRP conjugate. The assay sample and buffer are incubated together with TRY-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the TRY concentration since TRY from samples and TRY-HRP conjugate compete for the anti-TRY antibody binding site. Since the number of sites is limited, as more sites are occupied by TRY from the sample, fewer sites are left to bind TRY-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The TRY concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.