

Canine A Kinase Anchor Protein 1 ELISA Kit | AKAP1 elisa kit

Canine A Kinase Anchor Protein 1 ELISA Kit

CATALOG # AAA97245	REACTIVITY Canine
------------------------------	-----------------------------


PRODUCT SYNONYMS

A Kinase Anchor Protein 1; Canine A Kinase Anchor Protein 1 ELISA Kit; AKAP1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	2.5-50ng/mL
Sensitivity	0.1ng/mL
Intended Uses	This AKAP1 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine AKAP1. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Product Categories/Family for AKAP1 elisa kit	Signal Transduction
---	---------------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	4502015
NCBI GeneID	8165
NCBI Accession #	NP_003479.1
NCBI GenBank Nucleotide #	NM_003488.3

Molecular Weight	62,921 Da
NCBI Official Full Name	A-kinase anchor protein 1, mitochondrial
NCBI Official Synonym Full Names	A kinase (PRKA) anchor protein 1
NCBI Official Symbol	AKAP1
NCBI Official Synonym Symbols	AKAP; PRKA1; AKAP84; TDRD17; AKAP121; AKAP149; D-AKAP1; PPP1R43; SAKAP84
NCBI Protein Information	A-kinase anchor protein 1, mitochondrial; AKAP 149; D-AKAP-1; protein kinase A1; tudor domain containing 17; A-kinase anchor protein 149 kDa; protein kinase A anchoring protein 1; spermatid A-kinase anchor protein 84; protein phosphatase 1, regulatory subunit 43; dual-specificity A-kinase anchoring protein 1
UniProt Protein Name	A-kinase anchor protein 1, mitochondrial
UniProt Gene Name	AKAP1
UniProt Synonym Gene Names	AKAP149; PRKA1; AKAP 149; D-AKAP-1; PRKA1; S-AKAP84
UniProt Entry Name	AKAP1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/akap1-canine-elisa-kits/a-kinase-anchor-protein-1/97245>