

## Canine Midregional proAtrial Natriuretic Peptide ELISA Kit | MR-proANP elisa kit

Canine Midregional proAtrial Natriuretic Peptide ELISA Kit

CATALOG #  
**AAA97254**

REACTIVITY  
**Canine**


### PRODUCT SYNONYMS

Midregional proAtrial Natriuretic Peptide; Canine Midregional proAtrial Natriuretic Peptide ELISA Kit; MR-proANP elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                                                                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity              | Canine                                                                                                                                                                                            |
| Samples                 | Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate                                                                                                                        |
| Sensitivity             | 1.0 pmol/L.                                                                                                                                                                                       |
| Intended Uses           | This MR-proANP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine MR-proANP. This ELISA kit for research use only, not for therapeutic applications! |
| Preparation and Storage | Store all reagents at 2-8 degree C.                                                                                                                                                               |

### IMAGES



Standard Curve (Sample)

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://stage.aaabiotech.com/mr-proanp-canine-elisa-kits/midregional-proatrial-natriuretic-peptide/97254>