

Monkey intestinal alkaline phosphatase ELISA Kit | CIAP elisa kit

Monkey intestinal alkaline phosphatase ELISA Kit

CATALOG #
AAA97351

REACTIVITY
Monkey


PRODUCT SYNONYMS

intestinal alkaline phosphatase; Monkey intestinal alkaline phosphatase ELISA Kit; CIAP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Monkey
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	0.1 ng/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CIAP elisa kit

Principle of the assay: CIAP ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-CIAP antibody and an CIAP-HRP conjugate. The assay sample and buffer are incubated together with CIAP-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the CIAP concentration since CIAP from samples and CIAP-HRP conjugate compete for the anti-CIAP antibody binding site. Since the number of sites is limited, as more sites are occupied by CIAP from the sample, fewer sites are left to bind CIAP-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The CIAP concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.