

Mouse Uveal autoantigen with coiled-coil domains and ankyrin repeats ELISA Kit | Uaca elisa kit

Mouse Uveal autoantigen with coiled-coil domains and ankyrin repeats ELISA Kit

CATALOG # AAA99583	REACTIVITY Mouse
------------------------------	----------------------------


PRODUCT SYNONYMS

Uveal autoantigen with coiled-coil domains and ankyrin repeats; Mouse Uveal autoantigen with coiled-coil domains and ankyrin repeats ELISA Kit; Uaca/Kiaa1561/Uveal autoantigen with coiled-coil domains and ankyrin repeats/Nuclear membrane-binding protein/Nucling; Uaca elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	1416
Assay Type	Competitive ELISA, Coated with Antigen
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	0.313-20ng/ml
Sensitivity	0.188ng/ml
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	59850762
NCBI GeneID	55075
NCBI Accession #	NP_060473.2
NCBI GenBank Nucleotide #	NM_018003.2
Molecular Weight	161,510 Da
NCBI Official Full Name	uveal autoantigen with coiled-coil domains and ankyrin repeats isoform 1
NCBI Official Synonym Full Names	uveal autoantigen with coiled-coil domains and ankyrin repeats
NCBI Official Symbol	UACA
NCBI Official Synonym Symbols	NUCLING
NCBI Protein Information	uveal autoantigen with coiled-coil domains and ankyrin repeats

UniProt Protein Name	Uveal autoantigen with coiled-coil domains and ankyrin repeats
UniProt Gene Name	UACA
UniProt Synonym Gene Names	KIAA1561
UniProt Entry Name	UACA_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/uaca-mouse-elisa-kits/uveal-autoantigen-with-coiled-coil-domains-and-ankyrin-repeats/99583>